Fox-Flow

Sheath Fluid

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REF	100-05000-01	5 L
REF	100-05000-02	10 L

1. Účel použití

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1. Intended Use

The sheath fluid Fox-Flow is intended for use in flow cytometers. The product is intended for professional laboratory use only.

2. Method Principle

To accurately measure optical properties of cells with a flow cytometer, cells have to pass through the uniformly bright centre of focused laser beams. Light collection optics is focused on the intersection point of cells with the laser beams to pick up fluorescence and scattered light from cells. One purpose of the fluidic system is to move the cells one by one precisely through the sensing region in a liquid stream in such a way that each cell is illuminated by the same amount of light from the lasers [1]. For precise positioning of cells in a liquid jet, the hydrodynamic focusing technique is used in most cytometers: The cells in suspension are injected by a thin tubing in a laminar flow of a sheath fluid. Thus, the sheath flow speeds up and the diameter of sheath and sample flow is decreased [2].

The sheath fluid Fox-Flow is developed to comply with the principle described above. The product is a filtered solution which exhibits low level of autofluorescence, thereby reducing potential clogging of the system and ensuring a high fluorescent signal to noise ratio.

3. Main Components

Sodium chloride, potassium chloride, disodium phosphate, monopotassium phosphate, 2-phenoxyethan-1-ol, sodium fluoride, water, ProClin™ 300. The product does not contain sodium azide.

The product is available in two variants:

Name	Volume	Catalogue Number
Fox-Flow 5L	5 L	100-05000-01
Fox-Flow 10L	10 L	100-05000-02

The Safety Data Sheet is available to download on website of the manufacturer.

4. Storage and Handling

Store at 2–25 °C. Use the solution at temperatures specified in an instrument manual. Expiry date of the unopened product is 24 months. After opening, the solution is stable for 60 days at 2–25°C. Do not freeze.

The batch number and expiry date are indicated on the packaging.

5. Procedure

The product is ready to use without any prior preparations and dilution.

6. References

- Cossarizza A. et al. Guidelines for the Use of Flow Cytometry and Cell Sorting in Immunological Studies (Second Edition). *Eur J Immunol* **49**:1457–1973 (2019).
- Crosland-Taylor P. A Device for Counting Small Particles suspended in a Fluid through a Tube. *Nature* **171**: 37–38 (1953).

7. Symbols

REF	Catalogue Number	IVD	In Vitro Diagnostic Medical Device
LOT	Batch Number	CE	CE Marking of Conformity
	Expiry Date	i	Consult Instruction for Use
	Temperature Limit		Manufacturer
	This Way Up	UDI	Unique Device Identifier
	Fragile		

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Version: 1
Date of Issue: 24.06.2024
Revision Date: -
www.ita-intertact.com