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**SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO.,
LTD.**

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【SPECIFIKACE BALENÍ】

50mL×1

【SPECIFIKACE POUŽITÍ】

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【REFERENCE】

1. de Paoli P, Reitano M, Battistin S, Castiglia C, Santini G. Enumeration of human lymphocyte subsets by monoclonal antibodies and flow cytometry: a comparative study using whole blood or mononuclear cells separated by density gradient centrifugation. *J Immunol Methods*. 1984;72:349-353.
2. Ashmore LM, Shopp GM, Edwards BS. Lymphocyte subset analysis by flow cytometry. Comparison of three different staining techniques and effects of blood storage. *J Immunol Methods*. 1989;118:209-215.
3. Renzi P, Ginns LC. Analysis of T cell subsets in normal adults: comparison of whole blood lysis technique to Ficoll-Hypaque separation by flow cytometry. *J Immunol Methods*. 1987;98:53-56.
4. Romeu MA, Mestre M, González L, et al. Lymphocyte immunophenotyping by flow cytometry in normal adults: comparison of fresh whole blood lysis technique, Ficoll-Paque separation and cryopreservation. *J Immunol Methods*. 1992;154:7-10.
5. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. *Clin Immunol Newslett*. 1990;10(4):49-55.
6. Kidd PG, Vogt RF, Jr. Report of the workshop on the evaluation of T-cell subsets during HIV infection and AIDS. *Clin Immunol Immunopathol*. 1989;52:3-9.
7. Landay AL, Muirhead KA. Procedural guidelines for performing immunophenotyping by flow cytometry. *Clin Immunol Immunopath*. 1989;52:48-60.
8. Yu Jing, Nie Liping, Gui Yaoting, Peng Liming. A Study on Erythrocyte Lysate of Flow Cytometer Whole Blood. *Medical Study Journal*. 2007;36(7):74-77.
9. Wang Jianzhong. A Study on Methodology of Blood Lymphocyte Phenotype Analysis Using Flow Cytometry *Chinese Journal of Laboratory Medicine*. 2000;23(4):203-206.
10. Xavier Bossuyt, Gerald E. Marti, and Thomas A. Fleisher. Comparative Analysis of Whole Blood Lysis Methods for Flow Cytometry. *Cytometry*. 1997;30:124-133.
11. P. Mene'ndez, O. Redondo, A. Rodriguez, M.C., et al. Comparison Between a Lyse-and-Then-Wash Method and a Lyse-Non-Wash Technique for the Enumeration of CD341 Hematopoietic Progenitor Cells. *Cytometry*. 1998;34:264-271.
12. *Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes; Tentative Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1992. NCCLS document H42-T.
13. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture: Approved Standard*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS document H3-A3.
14. Du Xiumin, Sun Xiaoming, Yin Geping, Qi Falian. Factors Influencing the Clinical Application of Flow Cytometry. *Clinical Laboratory Journal*. 2002;20(6):365-366.
15. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in healthcare settings. *MMWR*. 1988; 37:377-388.
16. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue: Tentative Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS

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【KONTAKTY】

Výrobce: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adresa: Mindray Building, Keji 12th Road South, High-tech industrial park, Nanshan, Shenzhen 518057, P.R.China

Web: www.mindray.com

E-mail: service@mindray.com

Tel: +86 755 81888998

Fax: +86 755 26582680

EC-Zástupce: Shanghai International Holding Corp. GmbH(Europe)

Adresa: Eiffestraße 80, 20537 Hamburg, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

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07- 2014


In Vitro diagnostic
medical device

Batch code

European
Conformity

Authorized representative in
the European Community

Use-by
date

Consult Instructions
for use

Temperature
limit

Manufacturer

Catalogue
number

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SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (dále jen Mindray) vlastní práva na duševní vlastnictví těchto Mindray produktů a tohoto návodu k použití. Tento návod může zmiňovat informace chráněné copyright nebo patentem a nezprostředkovává žádnou licenci na základě patentního práva nebo copyright Mindray nebo jiných firem. Mindray má v úmyslu zachovat informace obsažené v tomto manuálu jako důvěrné. Jakékoliv prozrazení informací z tohoto manuálu jakýmkoliv způsobem bez písemného souhlasu od Mindray je striktně zakázáno. Zveřejnění, změna, reprodukování, distribuce, pronájem, přizpůsobování, překládání nebo jakékoliv jiné práce odvozené od této příručky jsou přísně zakázány bez písemného souhlasu od Mindray.